

Features

- Accurate and user friendly
- simple to setup and easy to operate
- accept 10,20,30 and 50ml syringes
- Membrane keypad with Menu driven LCD



Specification

No. of Channels	One
Flow Rate	0.1 to 500ml/hr
Set Volume	0.1 to 999.9 ml
Interval Time	0 to 99.9 sec.
Accuracy	± 2 %
Input Supply	230V +/- 10%, 50Hz, Single phase AC, 50 Watt
Display	Single line 8 Character Liquid crystal Display
Motor	DC Stepper Motor (Cont. Duty)
Drive	PWM Based DC Stepper Motor Drive
Material of construction	Powder Coated Aluminium
Weight	3.5kgs

Contact US

RAVEL HITEKS PVT LTD.,
 150A, Electronics Industrial Estate, Perungudi, Chennai-
 600 096. India.
 Ph:044-2496 3241 / 3251 / 1004 / 0825,
 Fax : 91-44-4204 9599
 Email:sales@ravelhiteks.com
 Web: www.ravelhiteks.com